

Envisioned and Designed by World Class
Laser Researchers and Clinicians to achieve
3D cluster based Cellular Healing and Rejuvenation

Dual Wavelength
Dual core AI backed programs


Novo Knee

Novo Hair

Combo Mode



Novolase Array







Technical Specifications

Wavelength 810nm - Infra Red Laser (12 Diodes)

660nm - Red Laser (4 Diodes)

Power 2400 mW/ 400mW +/- 10%

Duration 5, 10 and 15 minutes

Mode Continous and Pulsed Output

Dimensions 190 x 135 x 55mm

Weight 1.2Kg

Power Supply 5V 6A DC Power Adapter

Features

Treatment Mode	Time	Wavelength
Novo Knee	5/10 Minutes	810nm
Novo Hair	5/10/15 Minutes	660nm
Combo Mode	5 Minutes	810nm+660nm

Clinical Applications



Musculoskeleral Pain & Sports Injuries

- Shoulder, Knee, Neck, Back & Hip Pain (Sciatica)
- Bursitis/Plantar Fascitis/Tendonitis
- Ankle Injury/ Elbow Injuries
- · ACL Injury/ Hamstring Injuries



Wound Healing

- Surgical Wound
- Diabetic Foot
- Non Healing Ulcerations
- Pressure Sores
- Bed Sores



Neurogenic Pain

- Fribromyalgia
- Trigeminal Neuralgia
- Post Herpetic Neuralgia



Oro Facial Pain & Cancer Care

- Oral Mucositis and Dysphagia
- Oral Submucous Fibrosis and Trismus
- Temporomandibular Joint Disorder
- Xerostomia / Salivary Rejuvenation



Lymphedema

- Post Mastectomy Lymphedema
- Lymphedema



Degenerative Conditions

- Arthritis
- Herniated Disc
- Spinal Stenosis

Quality Assurance & Safety Certifications













GOOD MANUFACTURING PRACTICE

Corporate Social Responsibility



Part of the money that you pay towards the purchase of any Novolase product, will go towards making an Innocent child's life better!

Novolase Technologies is a registered monthly donor with **UNICEF India**

Donor ID: 357190

24 x 7 Laser Helpline



Now order anything at the click of a button! From a Laser unit to a **Polyamide tip** or a **NovoGreen Dye** or a **Teeth Whitening Gel**!

Just scan this on whatsapp or use our 24 x 7 Chat bot on www.novolase.in

www.novolase.in